



Greenwood Emergency Vehicles

Greenwood Emergency Vehicles Stock  
 Chris Gagnon  
 530 John Dietsch Blvd  
 North Attleboro, MA 02763

Greenwood Emergency Vehicles  
 Chris Gagnon

**Exp. Date:** 09/09/2026  
**Quote No:** PINK-0002  
**Job/Order No:** 702423,702425,702426  
**CITIMEDIC:** WC7148 WC- Type 3 CitiMedic  
**INFO:** SPEC\_WC WC-Wheeled Coach Boilerplate  
**T7:** WC-7148 WC- Type 3 CitiMedic Gen2 148"  
**EMS EQUIP:** EMS-LOOSE Factory Ship Loose - EMS Equipment

02/16/2026 14:29:27

PART NO	S	DESCRIPTION	QTY	ID
		<b>== WC-Wheeled Coach Boilerplate - 1.026 04/05/22 ==</b>		<b>REV</b>
		<b>ORDER COORDINATOR</b>		<b>REV</b>
00-00-0023		Order Coordinator - Soupy Manivanh	1	REV
		<b>INFORMATION</b>		<b>REV</b>
00-01-1030		RFI - Request for Information	1	REV
		A request for information (RFI) is a standard business process whose purpose is to collect written information about the capabilities of various suppliers. Normally it follows a format that can be used for comparative purposes.		
		An RFI is primarily used to gather information to help make a decision on what steps to take next. RFIs are therefore seldom the final stage and are instead often used in combination with the following: request for proposal (RFP), request for tender (RFT), and request for quotation (RFQ). In addition to gathering basic information, an RFI is often used as a solicitation sent to a broad base of potential suppliers for the purpose of conditioning supplier's minds, developing strategy, building a database, and preparing for an RFP, RFT, or RFQ.		
00-01-4015		Other (Not Applicable)	1	REV
00-01-4117		Manufacturer Certified ISO 9001	1	REV
		<b>WARRANTY</b>		<b>REV</b>
00-02-8710		WC, Warranty, Conversion, 12 Month	1	REV
00-02-8720		WC, Warranty, Paint, 60 month Prorated, Standard	1	REV
00-02-8730		WC, Warranty, Structural, 20 Years, Std., Mods	1	REV
00-02-8740		WC, Warranty, Limited Electrical	1	REV
00-02-8750		WC, Warranty, Cabinet Construction	1	REV
		<b>== WC- Type 3 CitiMedic Gen2 148" - 1.026 04/05/22</b>		<b>REV</b>
		<b>==</b>		

<b>PART NO</b>	<b>S</b>	<b>DESCRIPTION</b>	<b>QTY</b>	<b>ID</b>
<b>CERTIFICATIONS REQUIRED</b>				
00-05-0050		This unit built in accordance with KKK-A-1822-F in effect on the date of order with exceptions as noted herein.	1	REV
00-05-0100		WC, This unit built in accordance with KKK-A-1822F CN 10 Cabinet Requirements This unit as specified meets all requirements of KKK-A1822-F Change notice 10.  All cabinets shall be labeled as to their capacity Rating.	1	REV
<b>01- CHASSIS REQUIREMENTS</b>				
<b>GM CHASSIS</b>				
01-01-1027		CM, 2026 Chevrolet, G3500, 6.6L, Gas Cutaway, 139" WB Estimated Auto throttle OEM Alternators OEM  NOTE: Due to changes made to the OEM chassis by the manufacturer -some options may have to be changed to fit the new configuration.	1	REV
<b>ALTERNATORS</b>				
01-03-1001		Alternator(s) - Standard OEM	1	REV
<b>SUSPENSION</b>				
01-07-3001		Rear Suspension, OEM Springs	1	REV
<b>HIGH IDLE</b>				
<b>FRONT END ALIGNMENT</b>				
01-17-7501		Front End Alignment, None- QC Check -Standard	1	REV
<b>OUTSIDE REAR VIEW MIRRORS</b>				
01-20-0104		Mirrors, OEM-STD	1	REV
<b>WHEEL COVERS</b>				
01-24-0211		CM, Wheel Covers, Phoenix, G3500, Installed with Air Max Valve Extenders, std	1	REV
<b>BASE CONVERSION</b>				
01-45-0436		CM, Conversion, Type 3, 148 Module, Chevrolet -23-1	1	REV
<b>CAB EXTENSION</b>				
02-01-2017		CM, 8" Cab Extension, Chevrolet, CitiMedic G3500 Cab Extension to utilize same flooring as Module floor. ECN # 2096	1	REV
02-01-2127		CM, Interior Headroom, 72 Inches	1	REV
02-01-4002		CM, Module, Type 3, 148" Body Length, 90" wide, Duraseam Doors, Cutaway 148" L x 90" W	1	REV
02-01-9570		CM, Aisle Space, 49", CityMedic Plus (48-49" actual)	1	REV

PART NO	S	DESCRIPTION	QTY	ID
<b>EXTERIOR MODULE CONSTRUCTION</b>				
03-03-1063		CM, Stone Guard, Front, Diamond Plate, Type 3, Standard Height	1	REV
<b>FUEL FILLS &amp; SPLASH PLATE</b>				
03-03-8043		CM, Housing Fuel Fill, Cast, Type III	1	REV
03-03-8311		CM, Fuel Fill Splash Plate, Stainless, 10in Wide - installed below fuel fill full perimeter to bottom of skirt rail to protect body- attach to body using double sided tape - no screws	1	REV
<b>REAR KICKPLATE</b>				
03-03-9072		CM, Rear Kick plate, Under-Ride Bumper, w/ Tag recess, No lights, CM+	1	REV
<b>SUB-FLOOR SYSTEM</b>				
03-04-0176		CM, Floor, Interior, Non Wood, Expanded PVC, Type 3	1	REV
<b>MODULE ENTRANCE</b>				
<b>MODULE DOORS, HANDLES &amp; HOLD OPENS</b>				
03-06-0118		CM, Duraseam Hinged Doors with Hidden Jambes - Magnetic door switches, Standard	1	REV
03-06-0721		CM, Hold Open Rear Entry Doors, (2) 6" Door Grabber - Installed so doors will open as wide as possible - Mount at Bottom of doors.	1	REV
03-06-0819		CM, Hold Open, Gas Strut, Curb Side Entry Door, 35#	1	REV
03-06-1072		CM, Patient Entry Door Handles, Trimark Black/Chrome Pull Handle w/"Push Button" With Round Push button Emergency release top and bottom of interior door panels.	1	REV
03-06-1078		CM, Rear Entry Door Trailing Latch, Side release Paddle	1	REV
03-06-4205		CM, Electrical Feed to Required Doors, Spring Protector (Per Door)	7	REV
<b>POWER DOOR LOCKS</b>				
03-06-6213		CM, Power Door Locks, All Module doors (except Batt) Wiring to Module Compartment and Access Doors (except Battery) for electric locks, wiring to CS wall for Switch.  Includes the Tie in to the Cab door lock switches.	1	REV
03-06-6216		CM, Switch, Stealth, for Power Door Locks In grille on passenger's side.	1	REV
03-06-6218		CM, Relay Control Circuit, Power Door Locks	1	REV
<b>ENTRY DOOR INNER PANELS</b>				
03-06-7343		CM, Entry Door Panels, Durasafe, S/S, w/ Chevrons & Red Reflective Logo - Durasafe Design - Aluminum Center Panels covered: white/red reflective tape - Logo to be cut in the lower entry door panel with a Red reflective BM material (to match the chevron on the lock box) attached to an aluminum panel behind the logo.	1	REV

PART NO	S	DESCRIPTION	QTY	ID
		<b>SIDE DOOR STEP</b>		<b>REV</b>
03-08-2165		CM, Side Entry Step Well, Type III, Aluminum Diamond Plate	1	REV
03-08-3074		CM, Light, Side Entry Step well, LED, 2", Whelen #TOCACCCR	1	REV
		<b>WINDOWS</b>		<b>REV</b>
03-09-4200		CM, Window, Upper, CS/Slider ,Rr/Fixed, Privacy Tint, PAN Curbside window with lower slider section Rear windows fixed Privacy Tint	1	REV
		<b>REAR BUMPER AND REAR STEP CONSTRUCTION</b>		<b>REV</b>
03-12-2601		CM, Rear Bumper, Chevy Cutaway, w/ holes for 2" LED Lights with Gator Grip on flip up step. holes in sides of bumper for 2" marker light	1	REV
03-12-2611		CM, BLACK Polyurethane Coating, Rear Bumper Supports	1	REV
		<b>INSULATION</b>		<b>REV</b>
03-13-3111		CM, Insulation, Spray Foam, Walls, Ceiling, Under Floor and Inside all Doors This shall be done at Outside service location	1	REV
		<b>RUB RAIL AND FENDER RINGS</b>		<b>REV</b>
03-15-6735		CM, Rub Rails, Skirt Line, Extruded "C" channel (E-one Style), 148"	1	REV
03-15-6741		CM, Tape, Reflective, White, For Extruded Skirt rails	1	REV
03-15-8372		CM, Fender Flare, Bright Finish	1	REV
		<b>SPLASH GUARDS AND RUNNING BOARDS</b>		<b>REV</b>
03-16-1728		CM, Mud Flaps, Rear, w/ Wheeled Coach Logo	1	REV
03-16-3725		CM, Running Boards, Chevy Cutaway, w/Gator Grip,	1	REV
		<b>CAB TO MODULE ACCESS</b>		<b>REV</b>
03-19-2385		CM, Bulkhead Cabinet, Type 3, walk-thru w/8" Cab extension - with 8" extended cab area - With walk thru door	1	REV
03-19-2887		CM, Bulkhead, Door, Sliding w/Sliding Window with sliding Lexan window	1	REV
		<b>LICENSE PLATE HOLDER</b>		<b>REV</b>
03-20-4081		CM, License Plate Holder, Cast LP0002, for Under ride Bumper in rear kick plate, centered License holder comes with (2) LED lights.	1	REV
		<b>EXTERIOR COMPARTMENT CONSTRUCTION</b>		<b>REV</b>
04-01-0510		CM, Compartments, Std, Floor Drop Down Compartment #2 - 2.5" drop Compartment #4 - 4" drop Compartment #6 (5.5" for spare tire)	3	REV

PART NO	S	DESCRIPTION	QTY	ID
04-01-0512		CM, Compartments, Sweep out Compartments #1 & #8 are standard sweep out	2	REV
04-01-3020		CM, Door Sill Protector, Stainless, All Comps	1	REV
04-01-5202		CM, Compartments, Diamond Plate Finish Per Comp Machine stamped louvers located above the floor line for venting.	4	REV
		Comp #2,4,6		
		Comp #1 to have Baffle vent		
04-01-5204		CM, Exterior Battery Compartment, Diamond Plate Finish <b>EXTERIOR COMPARTMENT DOORS</b>	1	REV REV
04-02-0761		CM, Compartment Comp Handle, Trimark, Oval Pull Handle, Black/Chrome <b>COMPARTMENT DOOR HOLD OPEN</b>	6	REV REV
04-02-7772		CM, Hold Open, Gas Strut, Ext Comp, 30lb. Compartments "#6", "#1", & "#2"	3	REV
04-02-7773		CM, Hold Open, Gas Strut, Comp Door Only, 110 Degree, 30lb, Comp #4 to allow easier access to the Stair chair.	1	REV
04-02-7774		CM, Hold Open, Gas Strut, RF Comp, 30lb  #8	1	REV
		<b>COMPARTMENT LIGHTING</b>		REV
04-03-1479		WC, Oxygen Light, 3" Surface mount Whelen LED, Clear Mounted in Oxygen Compartment Wall #3	1	REV
04-03-1503		CM, Compartment Light, 4" Optronic,#BUL43-CBK, LED, Each Compartment  1 - Access Panel Wall 1 2 - Wall #1 4 - Ceiling 6 - Wall #3 8- Wall #3 in access panel	5	REV
		<b>COMPARTMENT #1 - STREETSIDE FWD</b>		REV
04-06-0453		CM, Compartment #1, 3/4 High, CitiMedic.	1	REV
		<b>COMPARTMENT #2 - STREETSIDE FWD WHEELWELL</b>		REV
04-07-2136		CM, Compartment #2, Citimedic , Type 3	1	REV
04-07-2645		CM, Compartment #2, Shelf, ADP, Fixed for Electrical storage Screw in shelf that extends from wall 1 to the recessed suction and wall 2 to flush with the face of the suction dogleg	1	REV

PART NO	S	DESCRIPTION	QTY	ID
04-07-2647		CM, Compartment #2, Shelf, Adjustable,w/Alum Track, First Shelf	1	REV
04-07-9200		Add Dogleg in Comp #2 for Recessed Suction in Face of A/A	1	REV
<b>COMPARTMENT #4 - STREETSIDE AFT</b>				<b>REV</b>
04-09-3251		CM, Compartment #4, Full Height, Dual Depth, Citimedic - Full height backboard storage on the #4 compartment door. - The upper section is 8" D due to dogleg for interior cabinetry. Bottom section is 15" deep.	1	REV
		Floor drop of 4" for the entire compartment		
04-09-4992		CM, Compartment #4, Strap, Seat Belt Style, Each (1) on comp #4 door for Backboards  (1) from wall #1 to Wall #3 about 6" from wall 2 for the Stair chair.	2	REV
04-09-8016		CM, Compartment #4, Bracket, Backboard Mount on Door Diamond plate "C-style boxes mounted to the door of Comp #4  one at top and bottom to allow backboards or scoop stretcher to be slid in from the outboard side.	1	REV
<b>COMPARTMENT #6 - CURBSIDE AFT</b>				<b>REV</b>
04-11-1260		CM, Compartment #6, Spare Tire, With 5.5" Drop pan for spare tire.	1	REV
04-11-4544		CM, Compartment #6, Shelf, Custom, Adjustable, First Shelf	1	REV
<b>COMPARTMENT #8 - CURBSIDE UPPER FORWARD</b>				<b>REV</b>
04-14-2330		CM, Compartment #8/8.5, Split Doors, 8" Cab Extension, Upper 8 compartment, Lower 8 with inside /outside access.  Lower compartment for (2) battery	1	REV
04-14-8030		CM, Compartment #8, Shelf, Adjustable, First Shelf RF Comp "8"/cabinet "U8" area - Centered top to bottom. C8 inside access only-Drawer	1	REV
04-14-8035		CM, Compartment #8, Shelf, Additional, Adjustable, Each RF Comp "8"/cabinet "L8" area centered top to bottom.	1	REV
<b>INTERIOR TRIM AND FEATURES</b>				<b>REV</b>
05-01-0125		CM, Interior Cabinets, Duralite Construction, Thermolite, CN10 Compliant Circuit board mounting board will be made of PVC	1	REV
05-02-0222		CM, Cabinet Shelf, Interior, Wood, Additional, Adjustable Cabinet L3.5 and U2	2	REV

<b>PART NO</b>	<b>S</b>	<b>DESCRIPTION</b>	<b>QTY</b>	<b>ID</b>
05-02-0225		CM, Shelf, Interior, Non-Wood, Adjustable, in Cabinet:	4	REV
		(1) U4		
		(1) U3.5		
		(1) U6		
		(1) C6		
05-02-0227		CM, Shelf, Interior, Non-Wood, Adjustable, Additional Cabinet U4 (2) total	1	REV
		<b>HEADLINER</b>		<b>REV</b>
05-06-2870		CM, Headliner, Flat, PVC, Hinged trough cover, Std	1	REV
		<b>FLOORING</b>		<b>REV</b>
05-07-6005		CM, Flooring, Lonplate II, Gunmetal Gray #424	1	REV
05-07-6030		CM, Flooring, Roll up, 3-4" Approx....., No Toe Kick The roll up is not recessed it is installed on the cabinet wall and Squad bench.	1	REV
05-07-9892		CM, Floor, Thresholds, .063 Aluminum, Painted w/Black Poly-urea coating	1	REV
		<b>GRAB RAILS</b>		<b>REV</b>
05-10-1537		CM, Entry Door Grab Handles (3), Custom "L" Shape, 1", Yellow  Yellow Powder Coat Finish Stainless Steel Grab Handles  ILOS	1	REV
		<b>Ceiling Grab Rail - Center</b>		<b>REV</b>
05-10-1939		CM, Ceiling Grab Rail, Overhead, 84", 1" Diameter Handicap Style, Yellow Yellow Powder Coat Finish	1	REV
		<b>IV FLUID HANGERS</b>		<b>REV</b>
05-11-5712		CM, IV Hangers, Cast, Fold Down Recessed Each Location: IV1CS IV2SS	2	REV
		<b>BASEWALL CABINET:</b>		<b>REV</b>
05-11-8029		CM, Cabinet, Base wall, 8" ext, CPR Seat, w/Telemetry Area, Recessed suction Cabinet L3.5	1	REV
05-11-8040		CM, Cabinet, Street side Upper, 148 Module, w/CPR Seat , 72" HR  Cabinet U2, U1 vented wood door	1	REV

PART NO	S	DESCRIPTION	QTY	ID
05-11-8044		CM, Street side rear, with side seat, 72" HR Cabinet U4 with (2) adjustable shelves Cabinet U3.5	1	REV
05-13-4204		CM, Tray, A/A, Std Aisle No Bio waste, Poured, Black	1	REV
05-13-4497		CM, Door, U1 Hinged Wood, Vented  hinged on right.	1	REV
05-13-4521		CM, Door, U2 Restocking Gray Polycarbonate Slider, CN 10 Compliant	1	REV
<b>Interior Street side #3 - CPR Seat</b>				<b>REV</b>
05-14-1520		CM, CPR side seat, Fixed Backrest	1	REV
05-14-1521		CM, CPR, Lid, Hinged Flip Up With Trimark Latch	1	REV
05-14-4313		CM, Door, U3.5, OHO Top Hinged Gray Polycarbonate (CWTG) CN 10 Compliant	1	REV
05-14-4333		CM, Door, U4, OHO Dual Hinged Gray Polycarbonate (CWDX) CN 10 Compliant	1	REV
05-14-4337		CM, Door, L3.5, Restocking (Top) Sliding Gray Polycarbonate CN 10 Co with two gas strut hold opens.	1	REV
05-14-4391		CM, CPR Seat, "U" Barrier Bar, Padded, 1.5" SS, Ash Gray Mounted at the front of the CPR side seat vertically in line with the Padded head cushion on the "L" cabinets.	1	REV
05-14-4547		CM, Telemetry Tray, Poured, Black	1	REV
<b>SQUAD BENCH AREA</b>				<b>REV</b>
05-16-0099		CM, Cabinet, Squad Bench, CM+, w/Tower at Head Cardiac Monitor Tower at head of bench with (1) drawer pulling out toward the Aisle and (1) Bio waste drawer underneath pulling out toward the Aisle.	1	REV
05-16-0104		CM, Lid, Squad Bench, Single, No Bio-Waste lift up	1	REV
05-16-0108		CM, Squad Bench, Handle, Trimark, STD	1	REV
05-16-0114		CM, Strut, Gas, 60 LB, Installed	1	REV
05-16-0117		CM, CPR, Handle, Trimark, STD	1	REV
05-16-0130		CM, SB Tower Tray, Poured, Gray	1	REV
05-16-0131		CM, SB Tower Tray, Poured, Black	1	REV
<b>SEAT BELTS/ SAFETY NETS</b>				<b>REV</b>
05-16-0223		CM, CPR Seat, Seatbelt, Assy, (4 Point), Per4Max, Black, (1) Each, CN-8 Above the lid	1	REV
05-16-0231		CM, Squad Bench, Seatbelt, Assy, (4 Point), Per4Max, Black, (2) Each Above the lid On squad bench	1	REV

PART NO	S	DESCRIPTION	QTY	ID
<b>CURBSIDE REAR</b>				
05-16-9980		CM, Cabinet, Curbside Rear, In Rear wall, U6,C6 Gray OHO doors, w/two shelves (1) adjustable shelf in cabinet U6 (1) adjustable shelf in cabinet C6	1	REV REV
<b>BIOHAZARD-</b>				
05-17-2514		CM, Glove Butler-Cabinet, (3),Above C/S door	1	REV REV
<b>CURBSIDE RIGHT STACK STORAGE #8</b>				
05-18-0220		CM, RF ALS, Two Opening Cabinet w/Ctr drawer, HVAC at Top	1	REV
05-18-5168		CM, Doors, RF, U8, Dual Gen III OHO, Polycarbonate Doors, Gray	1	REV
05-18-5232		CM, RF ALS, C8 Center Drawer, 6" high x Full width x Full depth Metal Drawer with Solid Face for maximum drawer space	1	REV
05-18-5364		CM, Doors, RF, L8, Dual Gen III OHO, Polycarbonate Doors, Gray	1	REV
<b>ATTENDANT SEAT</b>				
05-19-7302		CM, Attendant's Seat, EVS 1780 "Comfort" Hi-Back w/Per4Max belting	1	REV
05-19-7320		CM, Color, Ash Gray, Per4Max Belt, Black	1	REV
05-19-7517		CM, Attendant's Seat Base, EVS Swivel 8 position	1	REV
<b>INTERIOR COLORS</b>				
05-20-5254		CM, Mica, Gray Pampas #4168-01, Gloss	1	REV
<b>UPHOLSTERY</b>				
05-20-6008		CM, Upholstery, Seamless, Ash Gray	1	REV
05-20-6028		CM, Backrest Cushion, EVS, Ash Gray, w/ WC Logo, Each	3	REV
<b>CABINET KICK PANELS</b>				
05-21-5130		CM, Kick plate, Brushed SS, Face of S/B & Base wall, 148 CM	1	REV
<b>CABINET LATCHES- INTERIOR</b>				
05-21-8076		CM, Latch, M1, 3/4" Material, 2" round, Non Locking, Chrome-CN 10 Certified  (1)Sharps/Waste (1)C8 (1)SB Drawer (1)U1 (3) Circuit Board Access doors	7	REV
05-21-8082		CM, Latch, Southco M1, 2" S/S, (3/8" Door) Non-Locking, CN 10-Rated 10lb  Glove caddy	1	REV
<b>06 - ELECTRICAL GENERAL</b>				
06-01-0147		CM, 12V Electrical system, RelayPlex system Flash Pattern- Quad Flash 75 - Per KKK	1	REV REV
<b>CAMERAS</b>				
06-06-1446		CM, Camera, (2), Backup & Interior, Roscoe with 7" monitor Mount (1) camera exterior on rear of module, off center toward driver's side.(if there's a center light)	1	REV REV

PART NO	S	DESCRIPTION	QTY	ID
		Mount (1) STSC169I Dome camera interior above rear entry doors, off center toward driver's side.		
		Mount 7" STSK4730B monitor: on windshield in place of Mirror		
		Contains		
		STSC130B REAR-MOUNTED CAMERA		
		STSM244 MIRROR MONITOR (7")		
		STSH343 33' HARNESS		
		STSC169 INTERIOR DOME CAMERA		
		STSH343 33' HARNESS		
<b>COMMUNICATION</b>				<b>REV</b>
06-06-5547		CM, Coax Cable, RG58/U, Each, Terminate: From module roof Dome #4 to Action Area standard - with 6' pigtails	3	REV
		From module roof Dome #5,6 to Front Console standard - with 6' pigtails		
07-00-0615		CM, Light, in Circuit Board Area, for Electrical Troubleshooting - with integral switch; wired constant hot.	1	REV
<b>POWER SOURCES</b>				<b>REV</b>
07-01-0220		CM, Power Source, 12VDC, 20A, Ignition/Shoreline Hot 20amp 12 volt DC circuit routed to two locations, (1) pre-wire coil and tagged in action area and (1) pre-wire coil and tagged behind driver's seat.	1	REV
07-01-0221		CM, Power Source, 12 VDC, 15 Amp (+-), Std Model only - Ignition hot, terminating in cab console - with 6 foot tails (hot and ground).	1	REV
<b>VOLTMETER</b>				<b>REV</b>
07-02-1325		CM, Alarm, Low Voltage, Audio/Visual Light in cab console and Buzzer in cab.	1	REV
<b>BATTERY SYSTEM</b>				<b>REV</b>
07-04-8816		CM, Battery Switch, RelayPlex Electronic w/5 min Timer, Chevy T3 - Activated thru OEM ignition switch, timer function upon ignition "off". - Momentary rocker switch, on driver's side of cab console, to function as timer shutoff and also reactivate timer. - ONLY Module power is turned on/off. - Chassis power is NOT turned off.	1	REV
<b>12VDC OUTLETS</b>				<b>REV</b>

<b>PART NO</b>	<b>S</b>	<b>DESCRIPTION</b>	<b>QTY</b>	<b>ID</b>
07-10-1060		CM, Power Outlet, Power Point Style, (2), (1) Action Area (1) RF ALS Cabinet on standard 20-amp Battery hot circuit. - (1) outlet in action area - (1) in RF Cabinet "U8", wall #2, upper right	1	REV
07-10-1064		CM, Power Outlet, Power Point Style, Each State Location: Telemetry Area wall #2 near wall #1  L8 cabinet wall #2 near wall #3 up high  Above the tower on the squad bench  Constant hot	3	REV
07-10-1119		CM, Power Outlet, Kussmaul, USB Dual Port, 5VDC, 4.8 Amp, switch panel  Located in the cab switch panel per Electrical Engineering	1	REV
07-10-1120		CM, Power Outlet, Kussmaul, USB Dual Port, 5VDC, 4.8 Amp, switch panel, Addition Location:  Telemetry Area Action Area Above the tower on the squad bench  Battery HOT	3	REV
		<b>FRONT CONSOLE</b>		<b>REV</b>
08-01-3329		CM, Console, Drivers Switch & Radio, Center, Aluminum, Black, Chevy cutaway Black polyurethane coated .090 Aluminum 10.75 x 12.5 x 7" deep storage area at rear See below for Optional panels	1	REV
08-01-4666		WC, Flat Switch Panel	1	REV
		<b>ATTENDANT CONTROL PANEL</b>		<b>REV</b>
08-03-4528		CM, Action Area Switch Panel, Eaton Esport switches, Relocate CS Wall	1	REV
		<b>SIREN ELECTRONIC - CONTROL HEADS / AMPS</b>		<b>REV</b>
09-03-1944		CM, Siren Electronic, Whelen 295HFSA7, California Title 13 compliant	1	REV
		<b>SIREN SPEAKERS</b>		<b>REV</b>
09-03-2754		CM, Speakers, (2), Whelen, SA-315, In Grille area, Chevy Cutaway	1	REV
09-03-3349		CM, Switch, Siren/Horn Thru Horn Ring	1	REV
		<b>BACK-UP ALARMS</b>		<b>REV</b>

PART NO	S	DESCRIPTION	QTY	ID
09-05-1134		CM, Back-Up Alarm, w/ Cutoff, Auto Reset	1	REV
09-06-0205		CM, Switch, Emergency Master One Emergency Master Switch to turn all warning lights on.	1	REV
09-40-5178		CM, Cab Headliner Removal Fee, Type 1/3 for Monitor installation	1	REV
		<b>WARNING LIGHT FLASHERS</b>		<b>REV</b>
09-80-0164		CM, Automatic Secondary, Park/Neutral, w/Momentary Switch, Auto Reset When unit is placed into Park, Automatically defaults to cut off all "Primary Only" Defined lighting, the circuit would automatically reset when the gear shift lever is placed back into gear.  There shall be an override switch to turn back on the "Primary Only" defined lighting.  This removes the Primary and Secondary Buttons or Switches.	1	REV
09-80-7003		CM, Circuit, Lights Alternating Pattern Alternating pattern. for a (7) light grouping across the front  Passenger corner light will be Phase 1 Outboard light on Passenger side will be Phase 2 Inboard light on Passenger side will be Phase 1 Center light will be phase 2 Primary only Inboard light on Drivers side will be Phase 1 Outboard light on Drivers side will be Phase 2 Drivers side corner light will be phase 1  Passenger upper grille light will be Phase 2 Primary Drivers upper grille light will be Phase 1 Primary  Passenger lower grille light will be Phase 2 Primary Drivers lower grille light will be Phase 1 Primary  Curbside Front intersection will be phase 2 Primary only Front upper corner will be Phase 1 Rear intersection will be Phase 1 Primary only Rear upper corner will be Phase 2  Street side Front intersection will be phase 1 Primary only Front upper corner will be Phase 2 Rear intersection will be Phase 2 Primary only Rear upper corner will be Phase 1	1	REV

PART NO	S	DESCRIPTION	QTY	ID
		Rear		
		Street side upper corner will be Phase 1		
		Center upper light shall be Phase 2		
		Curbside upper corner will be Phase 1		
		Street side thru window will be Phase 2		
		Curbside thru window will be Phase 2		
		All via RelayPlex		
		<b>WARNING LIGHTS</b>		
09-95-1421		CM, Light, Whelen 900 LED, Red, Clear Lens, 90RR5SCR, External Flasher	12	<b>REV</b> REV
		FRONT UPPER BODY LIGHTING		
		Front Body - Upper Curbside Outer on coolbar		
		Front Body - Upper Curbside Inner on coolbar		
		Front Body - Upper Street side Inner on coolbar		
		Front Body - Upper Street side Outer on coolbar		
		MODULE SIDES LIGHTING		
		Street side Body - Upper Forward		
		Street side Body - Upper Rear		
		Curbside Body - Upper Forward		
		Curbside Body - Upper Rear		
		MODULE REAR LIGHTING		
		Rear Body - Upper Street side Corner		
		Rear Body - Upper Curbside Corner		
		MODULE CENTER REAR WINDOW LIGHTING		
		Rear Body - Center Street side Window		
		Rear Body - Center Curbside Window		
		Pri/Sec		
09-95-1423		CM, Light, Whelen 900 LED, Amber, Clear Lens, 90AA5SCR, External Flasher	1	REV
		Rear-Upper Center		
		Pri/Sec		
09-95-1424		CM, Light, Whelen 900 LED, Clear, 90CC5SCR, External Flasher Front Center on Coolbar	1	REV
		Pri only		
09-95-1435		CM, Flange, Super 900 LED, #90FLANGC	14	REV

PART NO	S	DESCRIPTION	QTY	ID
		<b>LED Series - 700 (Flange Separate)</b>		<b>REV</b>
09-95-2382		CM, Light, Whelen 700 LED, Red, Clear Lens, 70R02SCR FRONT UPPER BODY LIGHTING (2) Front Body - Upper Curbside Outer corner (2) Front Body - Upper Street side Outer corner  Pri/Sec	8	REV
		FRONT FENDER INTERSECTION LIGHTING Front Fender Warning - Intersectors - LH Side Front Fender Warning - Intersectors - RH Side Pri only		
		REAR WHEEL WELL INTERSECTION LIGHTING Rear wheel well Warning - Intersectors - LH Side Rear wheel well Warning - Intersectors - RH Side  Pri only		
		<b>WHELEN FLANGES 7 SERIES LIGHTS (Flange Separate)</b>		<b>REV</b>
09-95-3040		CM, Flange, Whelen, Chrome, 700 Series Light (1) For 7 Series Lights-CHROME flange	6	REV
09-95-3042		CM, Housings, Whelen, 700 Intersection Lights, Pair, Chevy Cutaway #7ICV6LRC  Front fender	1	REV
		<b>WARNING LIGHTS -WHELEN LIN3 (Flange Separate)</b>		<b>REV</b>
09-95-4643		CM, Light, Whelen, LIN3, LED-Red, RSRO2ZCR w/Flg FRONT GRILL LIGHTING Grille - LH Upper Corner Grille - RH Upper Corner  LOWER FRONT GRILL LIGHTING Grille - LH Upper Corner Grille - RH Upper Corner  Pri only	4	REV
		<b>EXTERIOR AUTOMOTIVE LIGHTING</b>		<b>REV</b>
10-01-0202		CM, Taillight Pkg, Whelen, LED Stop/Tail, Turn & LED Backup, w/WC logo - Brake- Steady Burn - Turn Arrow – sequential arrow	1	REV

PART NO	S	DESCRIPTION	QTY	ID
		- 700 LED Backup Clear lens over all three lights		
<b>ICC/MARKER LIGHTS</b>				<b>REV</b>
10-01-2261		CM, Light, Red, 2", Marker, LED, Each in bumper	2	REV
10-01-2262		CM, Clearance Lights, Red, LED (3) On rear as clearance lights.	3	REV
10-01-2265		CM, Light, Clearance, Amber LED (3) On Front of module as clearance lights.	3	REV
10-01-2268		CM, Corner Cap Lts, Multi, LED, (2) Amber (2) Red w/Flashers - High intensity flashers;Primary/Secondary- As marker lights with headlights	1	REV
<b>FLOOD AND LOAD SYSTEMS</b>				<b>REV</b>
10-02-2103		CM, Light, LED, Whelen, Clear, 900EZ, 904SLC, (1) w/ Chrome Flange, Steady Burn - (2) each side as "Scene lights" - (2) on rear over rear doors as "Load Lights"	6	REV
10-02-2136		CM, Scene Lights, C/S, Activated with C/S Door Master Hot.	1	REV
10-02-2138		CM, Rear Side Scene Lts to Activate w/Rev Gear. Rear most scene lights to come on in reverse, Rear load lights are wired to reverse standard.	1	REV
<b>11 - INTERIOR LIGHTING</b>				<b>REV</b>
11-01-0206		CM, Dome Lights, LED, Whelen,White Flange,(3)S/S,(3)C/S, CM+ Three (3) street side in liner and three (3) curbside in liner equally spaced  Each side controlled by separate switch in action area standard. High/Off/Low  NOTE: NO VISABLE FASTENERS OR HOLES.	1	REV
11-01-0235		CM, Timer, Momentary Switch, 15 Minute, Constant Hot Mount switch on C/S wall at the head of the squad bench in the standard location. Switch to allow activation and also deactivation. Wire to street side dome lights in liner on high.	1	REV
11-01-2871		WC, No Lighting Recessed into Trough Cover	1	REV
<b>SPOTLIGHTS / HANDHELD LIGHTS</b>				<b>REV</b>
11-02-0104		CM, Handheld LED Spotlight, Novatech, Nighthawk 2000L-101, 100,000 CP	1	REV
<b>ATTENDANT LIGHT</b>				<b>REV</b>
11-03-2054		CM, Light, Kinequip #621224W, 12" LED, Each	1	REV

PART NO	S	DESCRIPTION	QTY	ID
<b>SHORELINE INLET</b>				
12-01-9204		CM, Super Auto Eject, 20A, White Cover above comp "#2"	1	REV REV
<b>OUTLETS - 125VAC</b>				
12-02-2200		CM, Outlets, 125 VAC, Duplex, (1) Action Wall, (1) Right Front Cabinet - (1) outlet in action area - (1) outlet in the RF Cabinet "U8", wall #2, upper right.	1	REV REV
12-02-2204		CM, Outlet, 125 VAC, Duplex, Additional, Each Location: Telemetry Area wall #2 near wall #1  L8 cabinet wall #2 near wall #3 up high  CS Wall Above the tower on the squad bench	3	REV
<b>BLOCKHEATERS</b>				
12-02-4001		Block Heater - With OEM Plug This option does not include a switch. (STANDARD).	1	REV REV
<b>INVERTER - 125VAC</b>				
12-03-1704		CM, Inverter/Charger, Vanner LifeSine Mounted on floor of cabinet U1.	1	REV REV
<b>AUXILIARY ENVIRONMENTAL SYSTEMS</b>				
13-02-0250		CM, Exhaust Fan, Single 12V Blower with "clamshell" vent on exterior rear of unit.	1	REV REV
13-02-4004		CM, Generation I, Steel Cool Bar cover	1	REV
<b>CITIMEDIC BRUSHLESS HVAC SYSTEM</b>				
13-02-4425		CM, Heat/AC, C3, Combo, Brushless , w/Coolbar with thermostat in action area.  12vdc External Condenser, in Coolbar mounted on front of module. Painted to match front of module.	1	REV REV
13-02-4800		CM, Thermostat, Master Control Digital CM5000, Brushless	1	REV
13-02-6023		Cover, Cool Bar, (5) 900 series lights, LED Lights	1	REV
13-02-6029		CM, Air Conditioning Face Plate, Brushed Stainless Steel with WC Logo	1	REV
<b>COT MOUNTING PROVISIONS</b>				
<b>Safety Hooks</b>				
14-01-1269		CM, Cot Mount, Stryker, Rear Safety Hook Installed	1	REV REV
<b>Cot Mounts</b>				

PART NO	S	DESCRIPTION	QTY	ID
14-01-3121		CM, Stryker Performance/Power Load, Add Floor Structure Plates & Wiring for Futu Additional floor structure for future installation of a Stryker Power Load or Performance Load w/ inductive charging System.	1	REV
14-01-3124		CM, Cut Floor and Install Stryker Floor Plate, 6390-700-001A, for Performance For Performance Load Cot mounts.  Note: Stryker Deutsch connection to be coiled up on the interior of the module.	1	REV
14-01-3138		CM, Power source for Induction Charging, 12V, 10 Amp constant hot	1	REV
14-01-4028		Center Mount	1	REV
<b>OXYGEN AND AIR SYSTEMS</b>				
14-02-2302		CM, Rack, Multiversal, Zico #QR-MV, W/ Straps (Uni-strut track mounting is not std)	1	REV
14-02-2304		CM, Track, Unistrut for Oxygen Rack	1	REV
14-02-3040		CM, Electric O2 System w/Manual Bypass & Amico Digital Monitor includes 12V solenoid, wiring, transducer, manual bypass valve, digital readout panel and 50 psi regulator	1	REV
14-02-4192		CM, O2 System, Ohio Style, (2 AA, 1 Vac Diss, 1 Liner, 1 C/S)  Ohio style outlets (2) oxygen outlets in action area (1) vacuum port in action area (1) oxygen outlet in liner at head of cot (1) oxygen outlet curbside wall at head of squad bench	1	REV
14-02-4930		CM, Lexan Flap, w/Hinge	1	REV
<b>SUCTION SYSTEM</b>				
14-03-0200		CM, Suction Pump, SSCOR 90206 mounted in compartment "2"	1	REV
14-03-0210		CM, Suction Holder, SSCOR 23002 with w/Canister Clip & Regulator Conforms to J3043 requirements	1	REV
14-03-3000		WC, Suction Outlet and Container located in Action area recess in recessed pocket	1	REV
<b>PAINT AND FINISH</b>				
15-00-0200		WC, Wheeled Coach Paint Process	1	REV
15-01-2001		RM & CM, Module Paint, Single Color, OEM White  Chevy Type 3 White: PPG Code - FBCH1002053	1	REV

PART NO	S	DESCRIPTION	QTY	ID
15-01-5000		Cab Paint - Single Color, OEM White	1	REV
		<b>PAINT STRIPING DESIGN</b>		REV
15-02-1200		Paint Belt, None	1	REV
		<b>GRAPHICS AND LETTERING</b>		REV
		<b>DECALS AND NOMENCLATURE</b>		REV
15-03-6440		CM, Decal, Roof Star, Installed 32" roof star of Life	1	REV
		<b>SHIP LOOSE ITEMS</b>		REV
15-05-0050		CM, Bracket, Spare Tire Mounting, Ship Loose	1	REV
74-09-0099		Indemnification Statement The purchaser agrees to defend, indemnify and hold Rev Ambulance Group harmless from any claims, costs (including actual attorneys' fees), damages and liabilities caused in whole or in part by any alteration or modification of, or changes or additions to the purchased products OR use of product for purposes it was not designed or intended for.	1	REV
		<b>== Factory Ship Loose - EMS Equipment - 1.026 04/05/22 ==</b>		REV